

Tomato + Minerals

Medium Detail Information Medium No.029MediumTomato + MineralsCompositionTomato 200g , grinded , boiled in 15-20 minutes and filtered before use

Glucose	10g/l
Yeast extract	10g/l
CaCO ₃	2g/l
Solution A	2.5ml/l
Solution B	2.5ml/l
Distilled water	1L

pH 6.9 ± 0.2 at 25°C

Solution A:	KH ₂ PO ₄	10g/l
	K ₂ HPO ₄	10g/l
	Distilled water	1 L
Solution B:	MgSO ₄ .7 H ₂ O	40g/l
	NaCL	2g/l
	FeSO ₄	2g/l
	MnSO ₄	2g/l
	Distilled water	1 L Comment